

## Technical Report

### Buffered Mobile Phases in Reversed-Phase Liquid Chromatography

Aqueous mobile phase preparation is perhaps the most critical factor in reversed-phase chromatography (RPC) method development of ionic analytes. This includes consideration of the effects of pH on analyte retention, what buffer to use, what concentration of buffer to use, the solubility of the buffer in the organic modifier, the effect of the buffer on detection, and probably other considerations as well.

#### Retention

Mobile phase pH fundamentally impacts retention of ionic analytes in RPC. In consideration of the ionic functional groups, the neutral (uncharged) form of the analyte is better retained. Thus, the dissociation properties of these ionic functional groups will affect the retention of the analyte. For acidic moieties (usu. carboxylates) a pH below the  $pK_a$  (within limits) enhances retention of the compound, while for basic moieties (usu. amines), a pH above the  $pK_a$  (within limits) enhances retention. Therefore retention of non-ionic analytes is effected minimally by mobile phase pH.

It is within a pH range near the  $pK_a$  of a given functional group that dramatic affects on retention is observed. This becomes apparent in consideration of a mathematical description of chemical dissociation.

For a weak acid:  $HA \leftrightarrow H^+ + A^-$

and 
$$K_a = \frac{[H^+][A^-]}{[HA]}$$

Therefore, by algebraic rearrangements including definition of  $pH = -\log [H^+]$ ,

$$pH = pK_a + \log \frac{[A^-]}{[HA]}$$

The corresponding equation for a weak base is

$$pH = pK_a + \log \frac{[B]}{[BH^+]}$$

These last two equations are commonly known as the Henderson-Hasselbach equation. An important point to note is that in the last term, the unprotonated species is in the numerator and the protonated species is in the denominator, regardless of whether consideration is made of an acid or base.

From the Henderson-Hasselbach, it can be seen that at a pH far removed from the  $pK_a$ , the ratio of unprotonated-to-protonated forms is disparate, and a small change in the pH has minimal affect on this ratio. Thus, at pHs far removed from the  $pK_a$  of the analyte, moderate change in pH won't significantly affect retention. However, at pHs near the  $pK_a$ , small changes in pH will produce significant changes in the ratio of the two species. Recall that it is the uncharged form of the analyte that is better retained. Therefore, changing the pH within a range of values sufficiently close to the  $pK_a$ , will dramatically affect retention. Usually this range is  $pK_a \pm 1.5$  pH units.

#### Buffer Selection

Choice of buffer is typically governed by what pH is desired. With reference again, to the above Henderson-Hasselbach equation, when the pH is near the  $pK_a$ , significant changes in the ratio of unprotonated-to-protonated species can take place without major changes in the pH. It is when  $pH = pK_a$ , that the system is able to accommodate changes in proton (or hydroxide) concentration with minimal changes in pH. Thus, the term "buffer": the pH is buffered. We thus use the term "buffer" in common speech of aqueous solutions, to actually mean pH-buffer.

Therefore in order to provide a pH-buffer, it is necessary to select an ionizable chemical species whose  $pK_a$  is close to the desired pH. If the  $pK_a$  is too far removed from the desired pH, no significant pH-buffering will occur, and reproducibility of retention times may suffer. It is generally recommended to use a buffer within a pH range of its  $pK_a \pm 0.5$  pH units. However, the  $pK_a \pm 1.0$  pH units is acceptable. Table 1 lists some common buffers used for small molecule HPLC and their useful pH-buffer range.

**Table 1. Common Buffers, their pH Ranges, and UV-Vis Range**

Buffer	$pK_a$ (25°C)	Maximum Buffer Range	UV Cutoff (nm)
TFA	0.3		210 (0.1%)
Phosphate, $pK_1$	2.1	1.1-3.1	< 200
Phosphate, $pK_2$	7.2	6.2-8.2	< 200
Phosphate, $pK_3$	12.3	11.3-13.3	< 200
Citrate, $pK_1$	3.1	2.1-4.1	230
Citrate, $pK_2$	4.7	3.7-5.7	230
Citrate, $pK_3$	6.4	5.4-7.4	230
Carbonate, $pK_1$	6.1	5.1-7.1	< 200
Carbonate, $pK_2$	10.3	9.3-11.3	< 200
Formate	3.8	2.8-4.8	210 (10mM)
Acetate	4.8	3.8-5.8	210 (10mM)
Ammonia	9.3	8.3-10.3	200 (10mM)
Borate	9.2	8.2-10.2	n/a
TEA	10.8	9.8-11.8	< 200

Sometimes a buffer species is used in a mobile phase, but the pH is well outside of its ability to buffer it. An example is ammonium acetate, pH ~7. In this case, solid ammonium acetate is dissolved in water with no further pH adjustment. Ammonium acetate does not buffer pH between pH 6-8. It is usually in such cases that the use of the ammonium acetate is to provide some pH control (but not buffer) while taking advantage of its volatility downstream. It is also not unusual to see TFA or formic acid used alone at pHs outside of their buffering ranges: TFA as an ion-pairing reagent, formate for reasons of its volatility. Again, there are other reasons for their selection including that they do provide for pH control given the concentration of the analyte as compared to the acid (TFA, formate) in the bulk mobile phase.

## Buffer Concentration

As indicated above, better buffering occurs when the pH is closer to the  $pK_a$ . Additionally, this buffer capacity is enhanced by a higher concentration of the buffer itself. Higher buffer concentrations will give more reproducible separations of compounds partially ionized at the pH of the mobile phase, by reducing local perturbations of the pH of the migrating analyte peak. Generally, a buffer concentration of 10-50mM is adequate for small molecules. Note that the buffer concentration will decline upon mixing with the organic modifier.

## Buffer Solubility

An immediate concern regarding buffer concentration is consideration of its solubility upon addition of the organic modifier. This is especially important when performing gradient separations. Solubilities in mixtures of aqueous solution and a given organic solvent are not readily predictable. A suggested approach is to determine it empirically: try a given volume fraction of the buffer of interest (at a given concentration) with a given volume fraction of the organic solvent. Thoroughly mix the two phases and observe whether the solution remains transparent. If there is a solubility problem, a precipitate will form, or the solution may turn opaque. This sort of test must be done for the final conditions of a gradient run, to avoid consequential instrument and column problems related to buffer precipitation. While this is not fool-proof, a general rule is no more than 50% organic should be used with a buffer. However, keep in mind, this will depend on the specific buffer as well as its concentration.

## Buffer Effects on Detection

Often the choice of buffer is determined by the means of detection. For traditional UV detection, the buffer needs to be effectively transparent in this region. This too becomes more critical for gradient separations. All the buffers listed in Table 1 have low

enough absorption below 220nm, with the exception of citrate, which generally must be used above 230nm. Several of them can be used considerably lower than 220nm.

More common today, are issues related to compatibility with mass spectral (MS) detection. Though some common bench-top MS detectors have considerable rugged inlet systems, buffer volatility permits less intensive maintenance and thus less downtime of the MS detector. Preferred buffers addressing this issue of volatility are formate, acetate, carbonate, and ammonia. Furthermore, the issue of suppression of ionization is the more fundamental concern as regards to buffer selection in LC-MS method development. In this regard, formate and acetate are ideal choices for positive-ion mode detection. TFA, however, can negatively impact detector response even in positive-ion mode (see references 3 & 4), while with negative ion mode, strongly suppresses ionization. Acetic acid is good for negative-ion mode. Also for issues of mass detection, high concentrations of buffer can simply swamp-out any signal one might have wished to detect. In general, the less buffer needed the better, and 5-10mM is a suggested starting point. Thus, application of LC-MS does further limit buffer selection and buffer concentration. This often compromises the available options to develop an optimal separation, but is usually offset by the power of the MS.

### Practical Notes

- 1) In preparing the buffer, dissolve the solid in ~95% of the final volume desired. After the solution is adjusted to the pH desired, bring the solution to volume.
- 2) Always calibrate the pH meter either at the final adjusted pH (the final adjusted pH must match the pH of a standard) or calibrate the pH meter at values above and below the final adjusted pH. This is the only way to reliably measure pH with a pH-meter. You cannot reliably measure a value outside the range for which you have calibrated an instrument.
- 3) For consistent results, pH adjustments to the aqueous solution should be made before addition of organic. True, the actual pH will change after addition of the organic, but there is no good method for reliable measurement of pH after such addition. This way, at least the laboratory practice is consistent.
- 4) All buffer solutions should be filtered before use as mobile phase in HPLC. This is to remove any particulates that may already be in the water or introduced by the solid buffer when dissolved. A hydrophilic 0.45µ filter is recommended. Filtration should be a last step before use as is or before mixing with organic.

### References

1. Segel, I.H., Biochemical Calculations, John Wiley & Sons, 1976, New York.
2. Snyder, L.R., Kirkland, J.J., Glajch, J.L., Practical HPLC Method Development, John Wiley & Sons, 1997, New York.
3. Temesi, D., Law, B., 1999, The Effect of LC Eluent Composition on MS Response Using Electrospray Ionization, LC-GC, 17:626.
4. Apffel, A. et. al. 1995. Enhanced Sensitivity for Peptide Mapping with Electrospray Liquid Chromatography-Mass Spectrometry in the Presence of Signal Suppression Due to Trifluoroacetic Acid-Containing Mobile Phases., J. Chrom. A. 712:177.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com/supelco](http://www.sigma-aldrich.com/supelco)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 **AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 **AUSTRIA** • Sigma-Aldrich Handels GmbH • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • B-2880 Bornem **BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • 01239-010 São Paulo, SP **CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • 186 00 Praha 8 **DENMARK** • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand **FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex **GERMANY** • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen **GREECE** • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • H-1067 Budapest **INDIA** • Sigma-Aldrich Co. • Bangalore 560 048 **IRELAND** • Sigma-Aldrich Ireland Ltd. • Dublin 24 **ISRAEL** • Sigma Israel Chemicals Ltd. • Rehovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • 20151 Milano **JAPAN** • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 **KOREA** • Sigma-Aldrich Korea • Seoul **MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 58200 Kuala Lumpur  
**MEXICO** • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca **NETHERLANDS** • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht **NORWAY** • Sigma-Aldrich Norway • Torshov • N-0401 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań **PORTUGAL** • Sigma-Aldrich Química, S.A. • Sintra 2710 **RUSSIA** • Sigma-Aldrich Russia • Moscow 103062 **SINGAPORE** • Sigma-Aldrich Pte. Ltd.  
**SOUTH AFRICA** • Supelco • CH-9471 Buchs **SPAIN** • Sigma-Aldrich Química, S.A. • 28100 Alcobendas, Madrid **SWEDEN** • Sigma-Aldrich Sweden AB • 135 70 Stockholm  
**SWITZERLAND** • Supelco • CH-9471 Buchs **UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: [supelco@sial.com](mailto:supelco@sial.com)

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.